

## Antibodies Product Datasheet

<b>Product Name:</b>	VWF Antibody
<b>Catalog Number:</b>	C41995
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	IHC-p, IF(paraffin section)
<b>Dilutions:</b>	IHC-p 1:50-300
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from part region of human protein
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	VWF F8VWF
<b>Protein Name:</b>	von Willebrand factor (vWF) [Cleaved into: von Willebrand antigen 2 (von Willebrand antigen II)]
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	P04275
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	VWF Polyclonal Antibody detects endogenous levels of protein.
<b>Storage and Stability:</b>	-20°C/1 year