

Antibodies Product Datasheet

Product Name:	HERC5 Antibody
Catalog Number:	C48826
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HERC5
Concentration:	1 mg/ml
Gene Name:	HERC5 CEB1 CEBP1
Protein Name:	E3 ISG15--protein ligase HERC5 (EC 6.3.2.-) (Cyclin-E-binding protein 1) (HECT domain and RCC1-like domain-containing protein 5)
Other Name:	
Human Swiss Prot No.:	Q9UII4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HERC5 at Human
Storage and Stability:	Store at -20° C/1 year