

Antibodies Product Datasheet

Product Name:	POLR2L Antibody
Catalog Number:	C50632
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human POLR2L
Concentration:	1 mg/ml
Gene Name:	POLR2L
Protein Name:	DNA-directed RNA polymerases I, II, and III subunit RPABC5 (RNA polymerases I, II, and III subunit ABC5) (DNA-directed RNA polymerase III subunit L) (RNA polymerase II 7.6 kDa subunit) (RPB7.6) (RPB10 homolog)
Other Name:	
Human Swiss Prot No.:	P62875
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of POLR2L at Human, Mouse
Storage and Stability:	Store at -20° C/1 year