

Antibodies Product Datasheet

Product Name:	PBEF Antibody
Catalog Number:	C30929
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PBEF. at AA range: 411-460
Concentration:	1 mg/ml
Gene Name:	NAMPT
Protein Name:	PBEF
Other Name:	Nicotinamide phosphoribosyltransferase (NAmPRTase) (Nampt) (EC 2.4.2.12) (Pre-B-cell colony-enhancing factor 1) (Pre-B cell-enhancing factor) (Visfatin)
Human Swiss Prot No.:	P43490
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PBEF Polyclonal Antibody detects endogenous levels of PBEF
Storage and Stability:	-20°C/1 year