

Antibodies Product Datasheet

Product Name:	GATA-2/3 Antibody
Catalog Number:	B0933
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GATA3. AA range:274-323
Concentration:	1 mg/ml
Gene Name:	GATA2/GATA3
Protein Name:	Trans-acting T-cell-specific transcription factor GATA-2/3
Other Name:	GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3; GATA2; Endothelial transcription factor GATA-2; GATA-binding protein 2
Human Swiss Prot No.:	P23771/P23769
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GATA-2/3 Polyclonal Antibody detects endogenous levels of GATA-2/3 protein.
Storage and Stability:	-20°C/1 year