

Antibodies Product Datasheet

Product Name:	PP5 mouse mAb
Catalog Number:	N0834
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human PP5 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ppp5c
Protein Name:	
Other Name:	FLJ36922;FLJ55954;OTTHUMP00000165834;PP-T;PP5;PPP5;PPP5_HUMAN;PPP5C;PPT;Protein phosphatase 5, catalytic subunit;Protein phosphatase T;Serine/threonine protein phosphatase 5;Serine/threonine-protein phosphatase 5.
Human Swiss Prot No.:	P53041
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PP5 and does not cross-react with related protein
Storage and Stability:	-20°C/1 year