
Antibodies Product Datasheet

Product Name:	OR56B4 Antibody
Catalog Number:	R12-3225
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	OR56B4 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human OR56B4.
Concentration:	1 mg/ml
Gene Name:	OR56B4
Other Name:	Olfactory receptor 56B4; OR56B4
Human Swiss Prot No.:	Q8NH76
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	OR56B4 Antibodyantibody detects endogenous levels of OR56B4.
Storage and Stability:	Stable at -20°C for at least 1 year.