
Antibodies Product Datasheet

Product Name:	OSCAR Antibody
Catalog Number:	C45353
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	OSCAR
Protein Name:	Osteoclast-associated immunoglobulin-like receptor (Osteoclast-associated receptor) (hOSCAR) (Polymeric immunoglobulin receptor 3) (PIgR-3) (PIgR3) (Poly-Ig receptor 3)
Other Name:	
Human Swiss Prot No.:	Q8IYS5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	OSCAR Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year