
Antibodies Product Datasheet

Product Name:	ApoA-V(ab) Monoclonal Antibody
Catalog Number:	N0034
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ApoA-V expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	APOA5
Protein Name:	Apolipoprotein A-V
Other Name:	APOA5; RAP3; Apolipoprotein A-V; Apo-AV; ApoA-V; Apolipoprotein A5; Regeneration-associated protein 3
Human Swiss Prot No.:	Q6Q788
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	ApoA-V(ab) Monoclonal Antibody detects endogenous levels of ApoA-V(ab) protein.
Storage and Stability:	-20°C/1 year