

## Antibodies Product Datasheet

<b>Product Name:</b>	IDH3A Antibody
<b>Catalog Number:</b>	C47455
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from the Internal region of human IDH3A.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	IDH3A
<b>Protein Name:</b>	Isocitrate dehydrogenase [NAD] subunit alpha mitochondrial
<b>Other Name:</b>	IDH3A; Isocitrate dehydrogenase [NAD] subunit alpha, mitochondrial; Isocitric dehydrogenase subunit alpha; NAD(+)-specific ICDH subunit alpha
<b>Human Swiss Prot No.:</b>	P50213
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	IDH3A Polyclonal Antibody detects endogenous levels of IDH3A protein.
<b>Storage and Stability:</b>	-20°C/1 year