

Antibodies Product Datasheet

Product Name:	NCX1 Antibody
Catalog Number:	C18827
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from NCX1 . at AA range: 270-350
Concentration:	1 mg/ml
Gene Name:	SLC8A1
Protein Name:	Sodium/calcium exchanger 1
Other Name:	SLC8A1; CNC; NCX1; Sodium/calcium exchanger 1; Na(+)/Ca(2+)-exchange protein 1
Human Swiss Prot No.:	P32418
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NCX1 Polyclonal Antibody detects endogenous levels of NCX1 protein.
Storage and Stability:	-20°C/1 year