

Antibodies Product Datasheet

Product Name:	GRAP2 mouse mAb
Catalog Number:	N1018
Reactivity:	Human
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant protein of human GRAP2.
Concentration:	1 mg/ml
Gene Name:	grap2
Protein Name:	
Other Name:	Adapter protein GRID; GADS; GADS protein; GRAP-2; Grap2; GRAP2_HUMAN; GRB-2-like protein; GRB2-related adapter protein 2; GRB2-related adaptor protein 2; GRB2-related protein with insert domain; GRB2L; GRBLG; GRBX; Grf-40; Grf40 adapter protein; Grf40; GR
Human Swiss Prot No.:	O75791
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GRAP2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year