
Antibodies Product Datasheet

Product Name:	CD18 Monoclonal Antibody
Catalog Number:	N0106
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CD18 expressed in E. Coli.
Concentration:	
Gene Name:	ITGB2
Protein Name:	Integrin beta-2
Other Name:	ITGB2; CD18; MF17; Integrin beta-2; Cell surface adhesion glycoproteins LFA-1/CR3/p150; 95 subunit beta; Complement receptor C3 subunit beta; CD antigen CD18
Human Swiss Prot No.:	P05107
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	CD18 Monoclonal Antibody detects endogenous levels of CD18 protein.
Storage and Stability:	-20°C/1 year