

Antibodies Product Datasheet

Product Name:	PIAS 3 Antibody
Catalog Number:	C0363
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PIAS3. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	PIAS3
Protein Name:	E3 SUMO-protein ligase PIAS3
Other Name:	PIAS3; E3 SUMO-protein ligase PIAS3; Protein inhibitor of activated STAT protein 3
Human Swiss Prot No.:	Q9Y6X2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PIAS 3 Polyclonal Antibody detects endogenous levels of PIAS 3 protein.
Storage and Stability:	-20°C/1 year