

## Antibodies Product Datasheet

<b>Product Name:</b>	ITIH4 (80k, Cleaved-Arg689) Antibody
<b>Catalog Number:</b>	L0401
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Modification:</b>	Cleaved Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human ITIH4 (80k, Cleaved-Arg689)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ITIH4 IHRP ITIHL1 PK120 PRO1851
<b>Protein Name:</b>	ITIH4 (80k, Cleaved-Arg689)
<b>Other Name:</b>	Inter-alpha-trypsin inhibitor heavy chain H4 (ITI heavy chain H4) (ITI-HC4) (Inter-alpha-inhibitor heavy chain 4) (Inter-alpha-trypsin inhibitor family heavy chain-related protein) (IHRP) (Plasma kallikrein sensitive glycoprotein 120) (Gp120) (PK-120) [CI
<b>Human Swiss Prot No.:</b>	Q14624
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human ITIH4 (80k, Cleaved-Arg689)
<b>Storage and Stability:</b>	-20°C/1 year