

Antibodies Product Datasheet

Product Name:	Estrogen Receptor α Antibody
Catalog Number:	C47429
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB: 1:500-1:2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of Estrogen Receptor α
Concentration:	1 mg/ml
Gene Name:	ESR1
Protein Name:	Estrogen receptor (ER) (ER-alpha) (Estradiol receptor) (Nuclear receptor subfamily 3 group A member 1)
Other Name:	Estrogen receptor (ER) (ER-alpha) (Estradiol receptor) (Nuclear receptor subfamily 3 group A member 1); YM3483
Human Swiss Prot No.:	P03372
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Estrogen Receptor α protein
Storage and Stability:	-20°C/1 year