

Antibodies Product Datasheet

Product Name:	RHG44 Antibody
Catalog Number:	C45899
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 430-510
Concentration:	1 mg/ml
Gene Name:	ARHGAP44 KIAA0672 RICH2
Protein Name:	Rho GTPase-activating protein 44 (NPC-A-10) (Rho-type GTPase-activating protein RICH2) (RhoGAP interacting with CIP4 homologs protein 2) (RICH-2)
Other Name:	
Human Swiss Prot No.:	Q17R89
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	RHG44 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year