

Antibodies Product Datasheet

Product Name:	CD48 Antibody
Catalog Number:	C30357
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD48. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	CD48
Protein Name:	CD48 antigen
Other Name:	CD48; BCM1; BLAST1; CD48 antigen; B-lymphocyte activation marker BLAST-1; BCM1 surface antigen; Leukocyte antigen MEM-102; TCT.1; CD48
Human Swiss Prot No.:	P09326
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD48 Polyclonal Antibody detects endogenous levels of CD48 protein.
Storage and Stability:	-20°C/1 year