

Antibodies Product Datasheet

Product Name:	AGPAT3 Antibody
Catalog Number:	R12-2460
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	AGPAT3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGPAT3.
Concentration:	1 mg/ml
Gene Name:	AGPAT3
Other Name:	1-acyl-sn-glycerol-3-phosphate acyltransferase gamma 1-AGP acyltransferase 3; 1-AGPAT 3; Lysophosphatidic acid acyltransferase gamma; LPAAT-gamma; 1-acylglycerol-3-phosphate O-acyltransferase 3 Q9NRZ7
Human Swiss Prot No.:	
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	AGPAT3 Antibody antibody detects endogenous levels of AGPAT3.
Storage and Stability:	Stable at -20°C for at least 1 year.