

---

## Antibodies Product Datasheet

<b>Product Name:</b>	IRS-1 (Phospho-Tyr896) Antibody
<b>Catalog Number:</b>	A8220
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human IRS-1 around the phosphorylation site of Tyr896. AA range:862-911
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	IRS1
<b>Protein Name:</b>	Insulin receptor substrate 1
<b>Other Name:</b>	IRS1; Insulin receptor substrate 1; IRS-1
<b>Human Swiss Prot No.:</b>	P35568
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-IRS-1 (Y896) Polyclonal Antibody detects endogenous levels of IRS-1 protein only when phosphorylated at Y896.
<b>Storage and Stability:</b>	-20°C/1 year