
Antibodies Product Datasheet

Product Name:	MAP3K2 (Phospho-Ser520) Antibody
Catalog Number:	A8131
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MAP3K2 (Phospho-Ser520)
Concentration:	1 mg/ml
Gene Name:	MAP3K2 MAPKKK2 MEKK2
Protein Name:	MAP3K2 (Phospho-Ser520)
Other Name:	Mitogen-activated protein kinase kinase kinase 2 (EC 2.7.11.25) (MAPK/ERK kinase kinase 2) (MEK kinase 2) (MEKK 2)
Human Swiss Prot No.:	Q9Y2U5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.130% sodium azide.
Specificity:	This antibody detects endogenous levels of MAP3K2 (Phospho-Ser520) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year