

## Antibodies Product Datasheet

<b>Product Name:</b>	Alix Antibody
<b>Catalog Number:</b>	C12718
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000, ELISA(peptide)1:5000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human Alix
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PDCD6IP AIP1 ALIX KIAA1375
<b>Protein Name:</b>	Alix
<b>Other Name:</b>	Programmed cell death 6-interacting protein (PDCD6-interacting protein) (ALG-2-interacting protein 1) (Hp95)
<b>Human Swiss Prot No.:</b>	Q8WUM4
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of human Alix
<b>Storage and Stability:</b>	-20°C/1 year