

Antibodies Product Datasheet

Product Name:	Anti-Human Myelin Basic Protein(MBP) (ABT353) Mouse Monoclonal Antibody
Catalog Number:	V0407
Clone ID:	ABT353
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	MBP
Gene Name:	P02686
Other Name:	GDB, Golli MBP, Golli MBP; myelin basic protein, Hemopoietic MBP, HMBPR, HUGO, MBP, MBP_CAVPO, MBP_HUMAN, MGC99675, MLD, Myelin A1 protein, Myelin A1 Protein, basic, Myelin basic protein, Myelin Deficient, Myelin membrane encephalitogenic protein, OTTHUMP00000163776, OTTHUMP00000174387, OTTHUMP00000174388, SHI, Shiverer, SP
Isotype	IgG1, Kappa
Localization:	Brain/ Hippocampus
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Myelin Basic Protein.
Storage and Stability:	-20°C/1 year