

Antibodies Product Datasheet

Product Name:	AChR α 10 Antibody
Catalog Number:	C16982
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHRNA10. AA range:394-443
Concentration:	1 mg/ml
Gene Name:	CHRNA10
Protein Name:	Neuronal acetylcholine receptor subunit alpha-10
Other Name:	CHRNA10; NACHRA10; Neuronal acetylcholine receptor subunit alpha-10; Nicotinic acetylcholine receptor subunit alpha-10; NACHR alpha-10
Human Swiss Prot No.:	Q9GZZ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AChR α 10 Polyclonal Antibody detects endogenous levels of AChR α 10 protein.
Storage and Stability:	-20°C/1 year