

Antibodies Product Datasheet

Product Name:	USP22 Antibody
Catalog Number:	C48418
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human USP22
Concentration:	1 mg/ml
Gene Name:	USP22 KIAA1063 USP3L
Protein Name:	Ubiquitin carboxyl-terminal hydrolase 22 (EC 3.4.19.12) (Deubiquitinating enzyme 22) (Ubiquitin thioesterase 22) (Ubiquitin-specific-processing protease 22)
Other Name:	
Human Swiss Prot No.:	Q9UPT9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of USP22 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year