

Antibodies Product Datasheet

Product Name:	HXK III Antibody
Catalog Number:	C10708
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Hexokinase-3. AA range:811-860
Concentration:	1 mg/ml
Gene Name:	HK3
Protein Name:	Hexokinase-3
Other Name:	HK3; Hexokinase-3; Hexokinase type III; HK III
Human Swiss Prot No.:	P52790
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HXK III Polyclonal Antibody detects endogenous levels of HXK III protein.
Storage and Stability:	-20°C/1 year