

Antibodies Product Datasheet

Product Name:	P4HA2 Antibody
Catalog Number:	C48176
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human P4HA2
Concentration:	1 mg/ml
Gene Name:	P4HA2 UNQ290/PRO330
Protein Name:	Prolyl 4-hydroxylase subunit alpha-2 (4-PH alpha-2) (EC 1.14.11.2) (Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-2)
Other Name:	
Human Swiss Prot No.:	O15460
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of P4HA2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year