

Antibodies Product Datasheet

Product Name:	CD148 Antibody
Catalog Number:	C30482
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PTPRJ. AA range:861-910
Concentration:	1 mg/ml
Gene Name:	PTPRJ
Protein Name:	Receptor-type tyrosine-protein phosphatase eta
Other Name:	PTPRJ; DEP1; Receptor-type tyrosine-protein phosphatase eta; Protein-tyrosine phosphatase eta; R-PTP-eta; Density-enhanced phosphatase 1; DEP-1; HPTP eta; Protein-tyrosine phosphatase receptor type J; R-PTP-J; CD148
Human Swiss Prot No.:	Q12913
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD148 Polyclonal Antibody detects endogenous levels of CD148 protein.
Storage and Stability:	-20°C/1 year