

Antibodies Product Datasheet

Product Name:	Cytokeratin 14 Antibody
Catalog Number:	R12-2107
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytokeratin 14. AA range:102-151
Concentration:	1 mg/ml
Gene Name:	KRT14
Protein Name:	Keratin type I cytoskeletal 14
Other Name:	KRT14; Keratin; type I cytoskeletal 14; Cytokeratin-14; CK-14; Keratin-14; K14
Human Swiss Prot No.:	P02533
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cytokeratin 14 Polyclonal Antibody detects endogenous levels of Cytokeratin 14 protein.
Storage and Stability:	-20°C/1 year