

Antibodies Product Datasheet

Product Name:	CD25 Monoclonal Antibody(Q22)
Catalog Number:	N1372
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC: 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD25
Concentration:	
Gene Name:	IL2RA
Protein Name:	Interleukin-2 receptor subunit alpha
Other Name:	IL2RA; Interleukin-2 receptor subunit alpha; IL-2 receptor subunit alpha; IL-2-RA; IL-2R subunit alpha; IL2-RA; TAC antigen; p55; CD25
Human Swiss Prot No.:	P01589
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous human CD25 protein.
Storage and Stability:	-20°C/1 year