
Antibodies Product Datasheet

Product Name:	p53 Antibody
Catalog Number:	R12-3288
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	p53 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p53.
Concentration:	1 mg/ml
Gene Name:	TP53
Other Name:	Antigen NY-CO-13; Cellular tumor antigen p53; Phosphoprotein p53; TP53; tumor protein p53; Tumor suppressor p53
Human Swiss Prot No.:	P04637
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	p53 Antibody antibody detects endogenous levels of p53.
Storage and Stability:	Stable at -20°C for at least 1 year.