

Antibodies Product Datasheet

Product Name:	Anti-Human CD163 (ABT290) Mouse Monoclonal Antibody
Catalog Number:	V0221
Clone ID:	ABT290
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CD163
Gene Name:	Q86VB7
Other Name:	C163A_HUMAN, CD 163, CD163, CD163 antigen, CD163 molecule, Hemoglobin scavenger receptor, M130, M130 antigen precursor, Macrophage associated antigen, MM130, OTTHUMP00000238617, OTTHUMP00000238618, OTTHUMP00000238619, OTTHUMP00000238620, SCAR11, Scavenger receptor cysteine rich type 1 protein M130, sCD163, Soluble CD163
Isotype	IgG1, Kappa
Localization:	Tonsil
Positive Tissue Control:	Membranous, Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CD163 protein.
Storage and Stability:	-20°C/1 year