

Antibodies Product Datasheet

Product Name:	Collagen V α 2 (Cleaved-Leu1229) Antibody
Catalog Number:	L0387
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Collagen V α 2 (Cleaved-Leu1229)
Concentration:	1 mg/ml
Gene Name:	COL5A2
Protein Name:	Collagen V α 2 (Cleaved-Leu1229)
Other Name:	Collagen alpha-2(V) chain
Human Swiss Prot No.:	P05997
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Collagen V α 2 (Cleaved-Leu1229, protein was cleaved amino acid sequence between 1229-1230)
Storage and Stability:	-20°C/1 year