

## Antibodies Product Datasheet

<b>Product Name:</b>	SLC25A25 Antibody
<b>Catalog Number:</b>	C49379
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human SLC25A25
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SLC25A25 APC3 KIAA1896 MCSC3 SCAMC2 UNQ549/PRO1106
<b>Protein Name:</b>	Calcium-binding mitochondrial carrier protein SCaMC-2 (Mitochondrial ATP-Mg/Pi carrier protein 3) (Mitochondrial Ca(2+)-dependent solute carrier protein 3) (Small calcium-binding mitochondrial carrier protein 2) (Solute carrier family 25 member 25)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q6KCM7
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of SLC25A25 at Human, Mouse,Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year