

Antibodies Product Datasheet

Product Name:	APPL1 Antibody
Catalog Number:	C15314
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APPL1. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	APPL1
Protein Name:	DCC-interacting protein 13-alpha
Other Name:	APPL1; APPL; DIP13A; KIAA1428; DCC-interacting protein 13-alpha; Dip13-alpha; Adapter protein containing PH domain; PTB domain and leucine zipper motif 1
Human Swiss Prot No.:	Q9UKG1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	APPL1 Polyclonal Antibody detects endogenous levels of APPL1 protein.
Storage and Stability:	-20°C/1 year