
Antibodies Product Datasheet

Product Name:	cAMP Monoclonal Antibody
Catalog Number:	N0090
Reactivity:	Species independent
Applications:	ELISA
Dilutions:	ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	cAMP, conjugated to KLH.
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	cAMP Monoclonal Antibody detects cAMP protein.
Storage and Stability:	-20°C/1 year