

Antibodies Product Datasheet

Product Name:	St3Gal-III Antibody
Catalog Number:	R12-2368
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	St3Gal-III Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human St3Gal-III.
Concentration:	1 mg/ml
Gene Name:	ST3GAL3
Other Name:	CMP-N-acetylneuraminate-beta-1,4-galactoside alpha-2,3-sialyltransferase; Beta-galactoside alpha-2,3-sialyltransferase 3; Alpha 2,3-ST 3; Gal beta-1,3(4) GlcNAc alpha-2,3 sialyltransferase; N-acetyllactosaminide alpha-2,3-sialyltransferase; ST3Gal III; ST Q11203
Human Swiss Prot No.:	
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	St3Gal-III antibody detects endogenous levels of St3Gal-III.
Storage and Stability:	Stable at -20°C for at least 1 year.