

## Antibodies Product Datasheet

<b>Product Name:</b>	PDE6G (Phospho-Thr22) Antibody
<b>Catalog Number:</b>	A9601
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human PDE6G (Phospho Thr22)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PDE6G PDEG
<b>Protein Name:</b>	PDE6G (Phospho Thr22)
<b>Other Name:</b>	Retinal rod rhodopsin-sensitive cGMP 3',5'-cyclic phosphodiesterase subunit gamma (GMP-PDE gamma) (EC 3.1.4.35)
<b>Human Swiss Prot No.:</b>	P18545
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human, Mouse PDE6G (Phospho Thr22)
<b>Storage and Stability:</b>	-20°C/1 year