
Antibodies Product Datasheet

Product Name:	ELF5 Antibody
Catalog Number:	R12-2745
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	ELF5 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELF5.
Concentration:	1 mg/ml
Gene Name:	ELF5
Other Name:	ETS-related transcription factor Elf-5; E74-like factor 5; Epithelium-specific Ets transcription factor 2; ESE-2; Epithelium-restricted ESE-1-related Ets factor; ESE2
Human Swiss Prot No.:	Q9UKW6
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ELF5 Antibody antibody detects endogenous levels of ELF5.
Storage and Stability:	Stable at -20°C for at least 1 year.