

Antibodies Product Datasheet

Product Name:	c-Myc (Acetyl-Lys148) Antibody
Catalog Number:	D12256
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human c-Myc (Acetyl Lys148)
Concentration:	1 mg/ml
Gene Name:	MYC BHLHE39
Protein Name:	c-Myc (Acetyl Lys148)
Other Name:	Myc proto-oncogene protein (Class E basic helix-loop-helix protein 39) (bHLHe39) (Proto-oncogene c-Myc) (Transcription factor p64)
Human Swiss Prot No.:	P01106
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat c-Myc (Acetyl Lys148)
Storage and Stability:	-20°C/1 year