

Antibodies Product Datasheet

Product Name:	NFAT5 (Phospho-Ser1197) Antibody
Catalog Number:	A1220
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NFAT5 around the phosphorylation site of Ser1197. AA range:1171-1220
Concentration:	1 mg/ml
Gene Name:	NFAT5
Protein Name:	Nuclear factor of activated T-cells 5
Other Name:	NFAT5; KIAA0827; TONEBP; Nuclear factor of activated T-cells 5; NFAT5; T-cell transcription factor NFAT5; Tonicity-responsive enhancer-binding protein; TonE-binding protein; TonEBP
Human Swiss Prot No.:	O94916
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-NFAT5 (S1197) Polyclonal Antibody detects endogenous levels of NFAT5 protein only when phosphorylated at S1197.
Storage and Stability:	-20°C/1 year