

Antibodies Product Datasheet

Product Name:	Adrenergic Receptor alpha-2C Antibody
Catalog Number:	R12-2459
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Adrenergic Receptor alpha-2C Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Adrenergic Receptor alpha-2C.
Concentration:	1 mg/ml
Gene Name:	ADRA2C
Other Name:	Alpha-2C adrenergic receptor; Alpha-2C adrenoceptor; Alpha-2C adrenoreceptor; Alpha-2 adrenergic receptor subtype C4; ADRA2C; ADRA2L2; ADRA2RL2
Human Swiss Prot No.:	P18825
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Adrenergic Receptor alpha-2C Antibody antibody detects endogenous levels of Adrenergic Receptor alpha-2C.
Storage and Stability:	Stable at -20°C for at least 1 year.