

## Antibodies Product Datasheet

<b>Product Name:</b>	TNF-a mouse mAb
<b>Catalog Number:</b>	N0945
<b>Reactivity:</b>	Human
<b>Applications:</b>	ELISA
<b>Dilutions:</b>	ELISA 1:10000-20000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	tnf-a
<b>Protein Name:</b>	
<b>Other Name:</b>	APC1; APC1 protein; Cachectin; DIF; Differentiation inducing factor; Macrophage cytotoxic factor; MCF; Necrosin; TNF a; TNF alpha; TNF; TNF Macrophage Derived; TNF Monocyte Derived; TNF Superfamily Member 2; TNF superfamily, member 2; TNF, macrophage deri
<b>Human Swiss Prot No.:</b>	P01375
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects TNF-a.
<b>Storage and Stability:</b>	-20°C/1 year