

Antibodies Product Datasheet

Product Name:	KIT Antibody
Catalog Number:	R12-2953
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	KIT Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KIT.
Concentration:	1 mg/ml
Gene Name:	KIT
Other Name:	SCFR; Proto-oncogene tyrosine-protein kinase Kit; c-kit; CD117 antigen
Human Swiss Prot No.:	P10721
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	KIT Antibody antibody detects endogenous levels of KIT.
Storage and Stability:	Stable at -20°C for at least 1 year.