

Antibodies Product Datasheet

Product Name:	PGAM2 Antibody
Catalog Number:	C51160
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal PGAM2
Concentration:	1 mg/ml
Gene Name:	PGAM2 PGAMM
Protein Name:	Phosphoglycerate mutase 2 (EC 3.1.3.13) (EC 5.4.2.1) (EC 5.4.2.4) (BPG-dependent PGAM 2) (Muscle-specific phosphoglycerate mutase) (Phosphoglycerate mutase isozyme M) (PGAM-M)
Other Name:	Phosphoglycerate mutase 2 (EC 3.1.3.13) (EC 5.4.2.1) (EC 5.4.2.4) (BPG-dependent PGAM 2) (Muscle-specific phosphoglycerate mutase) (Phosphoglycerate mutase isozyme M) (PGAM-M)
Human Swiss Prot No.:	P15259
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PGAM2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year