

Antibodies Product Datasheet

Product Name:	Tyk 2 (Phospho-Tyr292) Antibody
Catalog Number:	P12-1064
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Tyk 2 (Phospho-Tyr292) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human Tyk 2 around the phosphorylation site of Tyr292.
Concentration:	1 mg/ml
Gene Name:	TYK2
Other Name:	TYK2; Non-receptor tyrosine-protein kinase TYK2
Human Swiss Prot No.:	P29597
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Tyk 2 (Phospho-Tyr292) antibody detects endogenous levels of Tyk 2 only when phosphorylated at Tyr292.
Storage and Stability:	Stable at -20°C for at least 1 year.