

Antibodies Product Datasheet

Product Name:	EphB1 Monoclonal Antibody
Catalog Number:	N0228
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EphB1 (aa19-133) expressed in E. Coli.
Concentration:	
Gene Name:	EPHB1
Protein Name:	Ephrin type-B receptor 1
Other Name:	EPHB1; ELK; EPHT2; HEK6; NET; Ephrin type-B receptor 1; ELK; EPH tyrosine kinase 2; EPH-like kinase 6; EK6; hEK6; Neuronally-expressed EPH-related tyrosine kinase; NET; Tyrosine-protein kinase receptor EPH-2
Human Swiss Prot No.:	P54762
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EphB1 Monoclonal Antibody detects endogenous levels of EphB1 protein.
Storage and Stability:	-20°C/1 year