
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | CS025 Antibody |
| Catalog Number: | C43301 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB |
| Dilutions: | WB 1:1000-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthesized peptide derived from human CS025 |
| Concentration: | 1 mg/ml |
| Gene Name: | C19orf25 |
| Protein Name: | UPF0449 protein C19orf25 |
| Other Name: | |
| Human Swiss Prot No.: | Q9UFG5 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of human CS025 |
| Storage and Stability: | -20°C/1 year |