

Antibodies Product Datasheet

Product Name:	SNRNP200 Antibody
Catalog Number:	C49080
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SNRNP200
Concentration:	1 mg/ml
Gene Name:	SNRNP200 ASCC3L1 HELIC2 KIAA0788
Protein Name:	U5 small nuclear ribonucleoprotein 200 kDa helicase (EC 3.6.4.13) (Activating signal cointegrator 1 complex subunit 3-like 1) (BRR2 homolog) (U5 snRNP-specific 200 kDa protein) (U5-200KD)
Other Name:	
Human Swiss Prot No.:	O75643
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SNRNP200 at Human
Storage and Stability:	Store at -20° C/1 year