

Antibodies Product Datasheet

Product Name:	SYK (Phospho-Tyr348) Antibody
Catalog Number:	A0091sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho derived from human SYK (Phospho-Tyr348)
Concentration:	1 mg/ml
Gene Name:	SYK
Protein Name:	Tyrosine-protein kinase SYK (EC 2.7.10.2) (Spleen tyrosine kinase) (p72-Syk)
Other Name:	Tyrosine-protein kinase SYK (EC 2.7.10.2) (Spleen tyrosine kinase) (p72-Syk)
Human Swiss Prot No.:	P43405
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This detects endogenous levels of SYK (Phospho-Tyr348)
Storage and Stability:	-20°C/1 year